

 Kathodeweg 5
 Molenwerf 19

 1627 LK - Hoorn
 1911 DB - Uitgeest

 T + 31 229 285010
 T + 31 251 320418

 E info@onbeperktinbeweging.nl
 W www.onbeperktinbeweging.nl

## Goliath 3 heavy duty wheelchair elevator

The Goliath wheelchair elevator is the most powerful elevator that is currently available on the European market. This beautifully styled elevator was designed to easily lift up the heaviest category of (platform) wheelchairs into a car with the aid of the supplied hoisting device or hoisting yoke.

Its maximum hoisting capacity of 200 kg makes it possible to even lift up the heaviest (platform) wheelchairs, including accessories. This is subject to the elevator being mounted in a car with sufficient space to place the (platform) wheelchair.

Exactly at the centre of gravity a hoisting point has been fitted to the (platform) wheelchair. As a result the (platform) wheelchair remains stable during lifting and turning.

The Goliath elevator can be turned into or outside of a car in a comfortable manner, either manually or by means of an electric rotating unit.

## Technical data:

Maximum hoisting capacity	:200 kg
Input current	:40A
Own weight	:25 kg
CE marking	







Onbeperkt in beweging