


S AFE U SER-FRIENDLY P URCHASED BY THE V.A. E LECTRIC R ELIABLE A FFORDABLE R ESALE VALUE RETAINED M AKES PARKING EASY



www.handicapsinc.com



The Superarm is the perfect Lift to fit your busy life style.

We all have lives that we need to tend to; careers, families to raise, errands to run and a million other things that require our attention. The last thing you need to be spending your energy on is fighting your mobility equipment. Which mobility equipment you choose can impact the quality of your life.

| Features | Benefits |
|---|---|
| No platform to fall from | Safest lift on the market* |
| Compact design | Doesn't block windows or doorway Leaves more room inside the van |
| Safely sets the user down on inclines | Adaptable to multiple terrains |
| No Platform in the doorway | Leaves doorway open |
| Minimal modification to vehicle | Protects resale value of your vehicle |
| Runs off vehicle battery | No extra battery needed |
| Lift is rated at 600 lbs | Will lift most chairs, electrics & scooters |
| Exceeds V.A. Standards | Funding available through some organizations |
| All Electric-Low Amperage | Less drain on your electrical system |
| No platform landing space required | Less parking space - only wheel- chair landing space needed |
| Little or no vehicle modifica- tion required | More money in your pocket |
| Designed by the disabled | Understanding your needs from experience |
| Low maintenance with no hydraulics | Little or no repair costs. Depend- able in all weather no leaks or freeze ups. |
| 3 year warranty (5 year limited) | 40+ years of service standing behind our products |
| Manual crank down on the motor | Back up feature, Never stranded if your vehicle electrical system fails |
| 4 Models available Vans-Mini- vans Motor homes- Basement Model Motorhomes | Fits your life style needs. Making life accessable in vans and most minivans, and motorhomes. |

*NHTSA's National Center for Statistics and Analysis Sept. 1997 report



Handicaps, Inc 4335 S. Santa Fe Drive Englewood, CO 80110



Designed by the disabled, this lift is simple, quiet, userfriendly and strong. The Superarm is the safest lift on the market*. It is reliable, low maintenance and all electric. No platform to block your doors. Less Space is required for parking. Affordable requiring minimal modification so that the vehicle resale value is retained. Purchased by the V.A. for those who qualify. Allows you the freedom to choose vehicle conversions you want for your life style and comfort rather than conversions to fit your mobility equipment. Each lift comes with mountain climber webbing straps, which attach to the user's wheelchair. The webbing will allow the chair to tip slightly back, giving head clearance as the user passes through the door.



(800) 782-4335 (303) 781-2062 Fax: (303) 761-6811